

For Research Use Only



FcZero-rAb® Biotin Anti-Human CD107a / LAMP1 Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98266-2

Basic Information

Catalog Number:
Biotin-FcA98266-2

Size:
100ug , 500 ug/ml

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
EG2611

GenBank Accession Number:
NM_005561.4

GeneID (NCBI):
3916

ENSEMBL Gene ID:
ENSG00000185896

UNIPROT ID:
P11279-1

Full Name:
lysosomal-associated membrane
protein 1

Calculated MW:
45 kDa

Purification Method:
Protein A purification

CloneNo.:
242116E8

Recommended Dilutions:
FC: 0.25 ug per 10⁶ cells in a 100 µl
suspension

Applications

Tested Applications:
FC

Species Specificity:
human

Positive Controls:

FC : Jurkat cells,

Background Information

LAMP1 (CD107a) is a heavily glycosylated membrane protein enriched in the lysosomal membrane. LAMP1 is extensively glycosylated with asparagine-linked oligosaccharides which protect it from intracellular proteolysis (PMID: 10521503). Although LAMP1 is expressed largely in the endosome-lysosomal membrane of cells, it is also found on the plasma membrane (PMID: 16168398). Elevated LAMP1 expression at the cell surface has also been detected during platelet and granulocytic cell activation, as well as in some tumor cells (PMID: 29085473). LAMP1 functions to provide selectins with carbohydrate ligands. This protein has also been shown to be a marker of degranulation on lymphocytes such as CD8+ and NK cells and may also play a role in tumor cell differentiation and metastasis (PMID: 18835598; 29085473; 9426697).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ Jurkat cells were surface stained with 0.25 ug Biotin Anti-Human CD107a / LAMP1 Rabbit RecAb (Biotin-FcA98266-2, Clone:242116E8) (red), or 5 ul FcZero-rAb® Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9)(blue) and Streptavidin-PE Conjugate (PE-PF00030). Cells were not fixed.