

For Research Use Only



FcZero-rAb® Biotin Anti-Human PSMA/GCPII Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98279

Basic Information

Catalog Number:

Biotin-FcA98279

Size:

100ug, 500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0216

GenBank Accession Number:

NM_004476.3

GeneID (NCBI):

2346

UNIPROT ID:

Q04609

Full Name:

folate hydrolase (prostate-specific membrane antigen) 1

Calculated MW:

84 kDa

Purification Method:

Protein A purification

CloneNo.:

242026E6

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

-

Applications

Tested Applications:

FC, ELISA

Species Specificity:

human

Positive Controls:

FC : LNCaP cells,

Background Information

PSMA (Prostate-specific membrane antigen) is also named as FOLH1, FOLH, NAALAD1, PSM and belongs to the peptidase M28 family. PSMA is a 100-120 kDa integral transmembrane glycoprotein, considered to be a highly specific marker of the prostate gland, and has successfully been used as a marker of circulating prostatic epithelial cells (PMID:10074909; 15680901). It is involved in conversion of the major neurotransmitter (NAAG) to NAA and free glutamate. It has 8 isoforms produced by alternative splicing.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ LNCaP cells (left) or HeLa cells (right) were surface stained with 0.25 ug Biotin Anti-Human PSMA/GCPII Rabbit RecAb (Biotin-FcA98279, Clone: 242026E6) (red) or 0.25ug Isotype Control (blue), and Streptavidin-PE Conjugate (PE-PF00030). Cells were not fixed.