

For Research Use Only

# CoraLite® Plus 405 Anti-Mouse CD41 (MWReg30)



Catalog Number: **CL405-65189**

## Basic Information

Catalog Number:

CL405-65189

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG1, kappa

GenBank Accession Number:

BC120493

GeneID (NCBI):

16399

UNIPROT ID:

Q9QUM0

Full Name:

integrin alpha 2b

Purification Method:

Affinity purification

CloneNo.:

MWReg30

Excitation/Emission maxima wavelengths:

399 nm / 422 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse peripheral blood platelets were surface stained with 0.5 ug Coralite® Plus 405 Anti-Mouse CD41 (CL405-65189, Clone: MWReg30) (red), or 0.5 ug Isotype Control. Cells were not fixed.



1X10<sup>6</sup> mouse peripheral blood platelets were surface stained with 0.5 ug Coralite® Plus 405 Anti-Mouse CD41 (CL405-65189, Clone: MWReg30) or 0.5 ug Isotype Control. Cells were not fixed.