

For Research Use Only

CoraLite® Plus 405 Anti-Human HLA-DR (L243)



Catalog Number: **CL405-65218**

Basic Information

Catalog Number:
CL405-65218

Size:
100tests , 5 µl/test

Source:
Mouse

Isotype:
IgG2a, kappa

GenBank Accession Number:
BC032350

GeneID (NCBI):
3122

Full Name:
major histocompatibility complex,
class II, DR alpha

Calculated MW:
254 aa, 29 kDa

Purification Method:
Protein A purification

CloneNo.:
L243

Excitation/Emission maxima
wavelengths:
399 nm / 422 nm

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information


Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 405 Anti-Human HLA-DR (CL405-65218, Clone:L243) or unstained. Cells were not fixed. Monocytes were gated.



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 405 Anti-Human HLA-DR (CL405-65218, Clone:L243) or unstained. Cells were not fixed. Monocytes were gated.