

For Research Use Only

CoraLite® Plus 405 Anti-Human CD71 (OKT9)



Catalog Number: CL405-65235

Basic Information

Catalog Number: CL405-65235	GenBank Accession Number: BC001188	Purification Method: Affinity purification
Size: 100tests , 5 µl/test	GeneID (NCBI): 7037	CloneNo.: OKT9
Source: Mouse	Full Name: transferrin receptor (p90, CD71)	Excitation/Emission maxima wavelengths: 399 nm / 422 nm
Isotype: IgG1, kappa	Calculated MW: 85 kDa	

Applications

Tested Applications:
FC
Species Specificity:
Human

Background Information


Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ PHA-treated (day 3) human PBMCs were surface stained with 5 ul CoraLite® Plus 405 Anti-Human CD71 (CL405-65235, Clone:OKT9) (red) or unstained (blue). The black dashed line indicates unactivated PBMCs surface stained with 5 ul CoraLite® Plus 405 Anti-Human CD71. Cells were not fixed. Lymphocytes were gated.