

For Research Use Only

# CoraLite® Plus 405 Anti-Rat CD43 (W3/13)



Catalog Number: **CL405-65237**

## Basic Information

Catalog Number:

CL405-65237

Size:

100ug , 500 µg/ml

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

GeneID (NCBI):

24796

UNIPROT ID:

P13838

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Affinity purification

CloneNo.:

W3/13

Excitation/Emission maxima  
wavelengths:

399 nm / 422 nm

## Applications

Tested Applications:

FC

Species Specificity:

Rat

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> Wistar rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 405 Anti-Rat CD43 (CL405-65237, Clone: W3/13) or 0.5 ug Control Antibody. Cells were not fixed.