

For Research Use Only

CoraLite® Plus 405 Anti-Mouse CD94 (18D3)



Catalog Number: **CL405-65240**

Basic Information

Catalog Number:

CL405-65240

Size:

100ug , 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC120854

GeneID (NCBI):

16643

UNIPROT ID:

O54707

Full Name:

killer cell lectin-like receptor,
subfamily D, member 1

Purification Method:

Affinity purification

CloneNo.:

18D3

Excitation/Emission maxima
wavelengths:

399 nm / 422 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 405 Anti-Mouse CD94 (CL405-65240, Clone: 18D3) or 0.5 ug Control Antibody. Cells were not fixed.