

For Research Use Only

# CoraLite® Plus 405 Anti-Mouse CD27 Rabbit Recombinant Antibody

Catalog Number: CL405-98160



## Basic Information

Catalog Number:

CL405-98160

Size:

100ug , 500 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001033126.2

GeneID (NCBI):

21940

UNIPROT ID:

P41272

Full Name:

CD27 antigen

Calculated MW:

28kDa

Purification Method:

Protein A purification

CloneNo.:

241378E11

Excitation/Emission maxima

wavelengths:

399 nm / 422 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Background Information

CD27 (also known as TNFRSF7) is a type I glycoprotein expressed on some B cells and the majority of T cells, and is a member of the tumor necrosis factor (TNF) receptor family. CD27 is required for generation and long-term maintenance of T cell immunity (PMID: 11062504). It is a receptor for CD70 (CD27L). Ligation of CD27 by CD70 induces strong ubiquitination of TRAF and the activation of both canonical and non-canonical NF-kappaB pathways, as well as the JNK pathway (PMID: 19426224). CD27 may also play a role in apoptosis through association with SIVA1.

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> mouse splenocytes were surface stained with APC Anti-Mouse CD3 (17A2) (APC-65077, Clone: 17A2) and 0.25ug CoraLite® Plus 405 Anti-Mouse CD27 Rabbit Recombinant Antibody (CL405-98160, Clone:241378E11) or CoraLite® Plus 405 Rabbit IgG Isotype Control RecAb (C405-98136, Clone: 240953C9). Cells were not fixed.



1x10<sup>6</sup> mouse splenocytes were surface stained with 0.25ug CoraLite® Plus 405 Anti-Mouse CD27Rabbit Recombinant Antibody (CL405-98160, Clone: 241378E11) (red) or CoraLite® Plus 405 Rabbit IgG Isotype Control RecAb (CL405-98136, Clone: 240953C9) (blue). Cells were not fixed.