

For Research Use Only

CoraLite®488-conjugated IκB Alpha Polyclonal antibody



Catalog Number: CL488-10268

Basic Information

Catalog Number: CL488-10268	GenBank Accession Number: BC002601	Purification Method: Antigen affinity purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 4792	Recommended Dilutions: IF 1:50-1:500
Source: Rabbit	Full Name: nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, alpha	Excitation/Emission maxima wavelengths: 488 nm/515 nm
Isotype: IgG	Calculated MW: 36 kDa	
Immunogen Catalog Number: AG0345		

Applications

Tested Applications: IF	Positive Controls: IF : HeLa cells,
Species Specificity: human	

Background Information

Nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, alpha (NFKBIA, synonyms: IKBA, MAD-3, NFKBI). NFKB1 or NFKB2 is bound to REL, RELA , or RELB to form the NFKB complex. The NFKB complex is inhibited by I-kappa-B proteins (NFKBIA or NFKBIB), which inactivate NF-kappa-B by trapping it in the cytoplasm. Phosphorylation of serine residues on the I-kappa-B proteins by kinases (IKBKA, or IKBKB) marks them for destruction via the ubiquitination pathway, thereby allowing activation of the NF-kappa-B complex. Activated NFKB complex translocates into the nucleus and binds DNA at kappa-B-binding motifs such as 5-prime GGGRNNYCC 3-prime or 5-prime HGGARNYYCC 3-prime.

Storage

Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CoraLite® 488 IκB Alpha antibody (CL488-10268) at dilution of 1:200.