

For Research Use Only

CoraLite® Plus 488-conjugated TIA1 Polyclonal antibody



Catalog Number: CL488-12133

Featured Product

Basic Information

| | | |
|---|---|---|
| Catalog Number: CL488-12133 | GenBank Accession Number: BC015944 | Purification Method: Antigen affinity purification |
| Size: 100ul , Concentration: 1000 µg/ml by Nanodrop; | GeneID (NCBI): 7072 | Recommended Dilutions: IF 1:50-1:500 |
| Source: Rabbit | Full Name: TIA1 cytotoxic granule-associated RNA binding protein | Excitation/Emission maxima wavelengths: 488 nm / 515 nm |
| Isotype: IgG | Calculated MW: 214 aa, 24 kDa, 43 kDa | |
| Immunogen Catalog Number: AG2778 | Observed MW: 40 kDa | |

Applications

| | |
|--|--|
| Tested Applications: IF | Positive Controls: IF : L02 cells, |
| Species Specificity: human, mouse, rat | |

Background Information

TIA1, also named as p40-TIA-1, is involved in alternative pre-RNA splicing and regulation of mRNA translation by binding to AU-rich elements (AREs) located in mRNA 3' untranslated regions (3' UTRs). It possesses nucleolytic activity against cytotoxic lymphocyte target cells. TIA1 may be involved in apoptosis. Two isoforms of this protein exist - 41kDa and 42kDa. one of these was a missense variant (P362L) in TIA1. Similar to the ALS-related disease proteins TDP-43, hnRNP A1, and FUS, TIA1 is an RNA-binding protein containing a prionlike LCD and assembles into membrane-less organelles, including SGs. Postmortem neuropathology of five TIA1 mutations carriers showed a consistent pathological signature with numerous round, hyaline, TAR DNA-binding protein 43 (TDP-43)-positive inclusions. TIA1 mutations significantly increased the propensity of TIA1 protein to undergo phase transition. In live cells, TIA1 mutations delayed stress granule (SG) disassembly and promoted the accumulation of non-dynamic SGs that harbored TDP-43. Moreover, TDP-43 in SGs became less mobile and insoluble.

Storage


Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed L02 cells using CoraLite® Plus 488 TIA1 antibody (CL488-12133) at dilution of 1:200.