

For Research Use Only

CoraLite® Plus 488-conjugated DAZL Polyclonal antibody



Catalog Number:CL488-12633

Basic Information

Catalog Number: CL488-12633	GenBank Accession Number: BC027595	Purification Method: Antigen affinity purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 1618	Recommended Dilutions: IF 1:50-1:500
Source: Rabbit	Full Name: deleted in azoospermia-like	Excitation/Emission maxima wavelengths: 488 nm / 515 nm
Isotype: IgG	Calculated MW: 295 aa, 33 kDa	
Immunogen Catalog Number: AG3326	Observed MW: 33 kDa	

Applications

Tested Applications: IF	Positive Controls: IF : mouse testis tissue,
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed mouse testis tissue using CoraLite® Plus 488 DAZL antibody (CL488-12633) at dilution of 1:100.



Immunofluorescent analysis of (4% PFA) fixed mouse testis tissue using CoraLite® Plus 488 DAZL antibody (CL488-12633) at dilution of 1:100.