

For Research Use Only

CoraLite® Plus 488-conjugated ASNS Polyclonal antibody



Catalog Number: CL488-14681

Featured Product

Basic Information

Catalog Number:

CL488-14681

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6393

GenBank Accession Number:

BC014621

GeneID (NCBI):

440

UNIPROT ID:

P08243

Full Name:

asparagine synthetase

Calculated MW:

64 kDa

Observed MW:

55-64 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF

Species Specificity:

human

Positive Controls:

IF : SH-SY5Y cells,

Background Information

Asparagine synthetase [glutamine-hydrolyzing](ASNS) is also named as TS11. The asparagine synthetase gene, encoding the enzyme that catalyzes the synthesis of asparagine and glutamate using glutamine and aspartate, is a gene for which transcription is highly regulated by the nutritional status of the cell (PMID:15385533). ASNS has some isoforms with the MW of 55-64 kDa.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed SH-SY5Y cells using CoraLite® Plus 488 ASNS antibody (CL488-14681) at dilution of 1:200.