

For Research Use Only

CoraLite® Plus 488-conjugated TP53INP1 Polyclonal antibody

Catalog Number:CL488-17872



Basic Information

Catalog Number: CL488-17872	GenBank Accession Number: BC074868	Purification Method: Antigen affinity purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 94241	Recommended Dilutions: IF/ICC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q96A56	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Isotype: IgG	Full Name: tumor protein p53 inducible nuclear protein 1	
Immunogen Catalog Number: AG12358	Calculated MW: 240 aa, 27 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : A431 cells,
Species Specificity: human, mouse, rat	

Background Information

9


Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed A431 cells using CoraLite® Plus 488 TP53INP1 antibody (CL488-17872) at dilution of 1:200.