

For Research Use Only

# CoraLite® Plus 488-conjugated ZW10 Polyclonal antibody

Catalog Number: CL488-24561

Featured Product



## Basic Information

### Catalog Number:

CL488-24561

### Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG19138

### GenBank Accession Number:

BC128264

### GeneID (NCBI):

9183

### UNIPROT ID:

O43264

### Full Name:

ZW10, kinetochore associated, homolog (Drosophila)

### Calculated MW:

779 aa, 89 kDa

### Observed MW:

88-90 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

IF/ICC 1:50-1:500

### Excitation/Emission maxima wavelengths:

493 nm / 522 nm

## Applications

### Tested Applications:

IF/ICC

### Species Specificity:

human

### Positive Controls:

IF/ICC : HeLa cells,

## Storage

### Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

### Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using Coralite® Plus 488 ZW10 antibody (CL488-24561) at dilution of 1:200.