

For Research Use Only

CoraLite® Plus 488-conjugated Cardiac Troponin T Polyclonal antibody



Catalog Number: **CL488-26592**

Basic Information

Catalog Number: CL488-26592	GenBank Accession Number: BC002653	Purification Method: Antigen affinity purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 7139	Recommended Dilutions: IF 1:50-1:500
Source: Rabbit	Full Name: troponin T type 2 (cardiac)	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Isotype: IgG	Calculated MW: 36 kDa	
Immunogen Catalog Number: AG24911	Observed MW: 34-40 kDa	

Applications

Tested Applications: IF	Positive Controls: IF : mouse heart tissue,
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.**

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed mouse heart tissue using CoralLite® Plus 488 Cardiac Troponin T antibody (CL488-26592) at dilution of 1:200.