

For Research Use Only

# CoraLite® Plus 488-conjugated Integrin Alpha V Polyclonal antibody



Catalog Number: CL488-27096

Featured Product

## Basic Information

<b>Catalog Number:</b> CL488-27096	<b>GenBank Accession Number:</b> BC136442	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 100ul , Concentration: 1000 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 3685	<b>Recommended Dilutions:</b> IF 1:50-1:500
<b>Source:</b> Rabbit	<b>Full Name:</b> integrin, alpha V (vitronectin receptor, alpha polypeptide, antigen CD51)	<b>Excitation/Emission maxima wavelengths:</b> 493 nm / 522 nm
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 1048 aa, 116 kDa	
<b>Immunogen Catalog Number:</b> AG25825	<b>Observed MW:</b> 130 kDa	

## Applications

<b>Tested Applications:</b> IF	<b>Positive Controls:</b> IF : MCF-7 cells,
<b>Species Specificity:</b> human	

## Background Information

### Storage


**Storage:**  
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using CoraLite® Plus 488 Integrin Alpha V antibody (CL488-27096) at dilution of 1:100.