

For Research Use Only

CoraLite® Plus 488-conjugated CACNA2D1 Polyclonal antibody



Catalog Number: CL488-27453

Basic Information

Catalog Number:

CL488-27453

Size:

100ul , Concentration: 1000 µg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26190

GenBank Accession Number:

NM_000722

GeneID (NCBI):

781

UNIPROT ID:

P54289

Full Name:

calcium channel, voltage-dependent,
alpha 2/delta subunit 1

Calculated MW:

125 kDa

Observed MW:

125 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF

Species Specificity:

mouse

Positive Controls:

IF : mouse eye tissue,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse eye tissue using CoraLite® Plus 488 CACNA2D1 antibody (CL488-27453) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse eye tissue using CoraLite® Plus 488 CACNA2D1 antibody (CL488-27453) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).