

For Research Use Only

CoraLite®488-conjugated Anti-Mouse CD16 / CD32 (2.4G2)



Catalog Number: CL488-65080

Basic Information

Catalog Number:

CL488-65080

Size:

100 tests, 5 µl/test

Source:

Rat

Isotype:

IgG2b

GenBank Accession Number:

BC038070

GeneID (NCBI):

14130

Full Name:

Fc receptor, IgG, low affinity IIb

Purification Method:

Affinity purification

CloneNo.:

2.4G2

Excitation/Emission maxima wavelengths:

488 nm/515 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

CD16 (FcγRIII) is a 50-70 kDa low affinity Fc receptor found on the surface of natural killer cells, neutrophil polymorphonuclear leukocytes, monocytes and macrophages. CD16 mediates antibody-dependent cellular cytotoxicity (ADCC) and other antibody-dependent responses, such as phagocytosis. CD32 (FcγRII) is a 40 kD transmembrane glycoprotein that binds to the Fc region of IgG with low affinity. CD32 is present on phagocytic cells such as macrophages and neutrophils, and is involved in the process of phagocytosis and clearing of immune complexes. The 2.4G2 antibody reacts with the extracellular domains of mouse CD16 (FcRIII) and CD32 (FcRII) (PMID: 90108; 8406898).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ C57BL/6 mouse splenocytes were surface stained with APC-Anti-Mouse CD3ε (APC-65060, Clone: 145-2C11) and 5 ul CoraLite®488-conjugated Anti-Mouse CD16 / CD32 (CL488-65080, Clone: 2.4G2). Cells were not fixed.



1X10⁶ C57BL/6 mouse splenocytes were surface stained with APC-Anti-Mouse CD3ε (APC-65060, Clone: 145-2C11) and CoraLite®488-conjugated Rat IgG2b isotype control antibody. Cells were not fixed.