

For Research Use Only

CoraLite® Plus 488 Anti-Mouse CSF1R/CD115 (AFS98)

Catalog Number: CL488-65095



Basic Information

Catalog Number:

CL488-65095

Size:

100ug , 500 µg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC036343

GeneID (NCBI):

12978

Full Name:

colony stimulating factor 1 receptor

Purification Method:

Affinity purification

CloneNo.:

AFS98

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



1X10⁶ Balb/c mouse peritoneal macrophages were surface stained with 0.5 ug CoraLite® Plus 488-conjugated Anti-Mouse CD115 (c-fms) (CL488-65095, Clone: AF598) (green) or CoraLite® Plus 488-conjugated Rat IgG2a isotype control (black). Cells were not fixed.