

For Research Use Only

# CoraLite® Plus 488 Anti-Mouse CD24 (M1/69)



Catalog Number: **CL488-65121**

## Basic Information

Catalog Number:

CL488-65121

Size:

100ug , 500 µg/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC075622

GeneID (NCBI):

12484

Full Name:

CD24a antigen

Purification Method:

Affinity purification

CloneNo.:

M1/69

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 488 Anti-Mouse CD24 (CL488-65121, Clone:M1/69) (red), or 0.5 ug Control Antibody. Cells were not fixed.