

For Research Use Only

CoraLite® Plus 488 Anti-Human CD7 (4H9)



Catalog Number: **CL488-65203**

Basic Information

Catalog Number: CL488-65203	GenBank Accession Number: BC009293	Purification Method: Protein A purification
Size: 100tests , 5 µl/test	GeneID (NCBI): 924	CloneNo.: 4H9
Source: Mouse	Full Name: CD7 molecule	
Isotype: IgG2a, kappa	Calculated MW: 240 aa, 25 kDa	

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface co-stained with APC Anti-Human CD4 and 5 ul CoraLite® Plus 488 Anti-Human CD7 (CL488-65203, Clone:4H9) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 488 Anti-Human CD7 (CL488-65203, Clone:4H9) (red) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.