

For Research Use Only

CoraLite® Plus 488 Anti-Mouse CD81 (EAT2)



Catalog Number: **CL488-65220**

Basic Information

Catalog Number: CL488-65220	GenBank Accession Number: BC011433	Purification Method: Protein G purification
Size: 100ug , 500 µg/ml	GeneID (NCBI): 12520	CloneNo.: EAT2
Source: Armenian Hamster	Full Name: CD81 antigen	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Isotype: IgG		

Applications

Tested Applications:
FC

Species Specificity:
mouse

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with 1 ug CoraLite® Plus 488 Anti-Mouse CD81 (CL488-65220, Clone:EAT2), or 1 ug Isotype Control. Cells were not fixed.