

For Research Use Only

# CoraLite® Plus 488 Anti-Rat CD45 (OX1)



Catalog Number: **CL488-65228**

## Basic Information

Catalog Number:

CL488-65228

Size:

100ug, 0.5 mg/ml

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

XM\_006249910

GeneID (NCBI):

24699

UNIPROT ID:

P04157

Full Name:

protein tyrosine phosphatase, receptor type, C

Calculated MW:

143 kDa

Purification Method:

Affinity purification

CloneNo.:

OX1

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

FC

Species Specificity:

rat

## Background Information

CD45, also known as protein tyrosine phosphatase, receptor type C, is a type I transmembrane protein expressed on the surface of all haematopoietic cells with the exception of erythrocytes and platelets (PMID: 3489673; 28615666). CD45 is a pan-haematopoietic cell marker and has been shown to be essential for T- and B-cell activation and signalling (PMID: 9429890; 16378097).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 488 Anti-Rat CD45 (CL488-65228, Clone:OX1) (red), or 0.5 ug Isotype Control. Cells were not fixed.