

For Research Use Only

CoraLite® Plus 488 Anti-Mouse CD94 (18D3)



Catalog Number: **CL488-65240**

Basic Information

Catalog Number:

CL488-65240

Size:

100ug , 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC120854

GeneID (NCBI):

16643

UNIPROT ID:

O54707

Full Name:

killer cell lectin-like receptor,
subfamily D, member 1

Purification Method:

N/A

CloneNo.:

18D3

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 488 Anti-Mouse CD94 (CL488-65240, Clone:18D3), or 0.5 ug Isotype Control. Cells were not fixed.