

For Research Use Only

CoraLite® Plus 488 Mouse IgG2a Isotype Control (MOPC-173)



Catalog Number: CL488-65244

Basic Information

Catalog Number:

CL488-65244

Size:

100ug, 0.2 mg/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Affinity purification

CloneNo.:

MOPC-173

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC

Species Specificity:

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface stained with 0.2 ug CoraLite® Plus 488 Anti-Mouse IgG2a Isotype Control (CL488-65244, Clone: MOPC-173) (red), CL488 Anti-Human CD3 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.