

For Research Use Only

# CoraLite® Plus 488 Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody

Catalog Number: **CL488-65571**



## Basic Information

Catalog Number:

CL488-65571

Size:

100ug , 500 µg/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

MOPC-173

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

FC

Species Specificity:

n/a

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.

## Selected Validation Data



$1 \times 10^6$  human PBMCs were surface stained with 0.2  $\mu$ g Coralite® Plus 488 Mouse IgG2a Isotype Control Recombinant Antibody (CL488-65571, Clone: MOPC-173) (red), or Coralite® Plus 488 Anti-Human CD8 (UCHT4) (CL488-65204, Clone: UCHT4) (green), or Coralite® Plus 488 Mouse IgG2a Isotype Control (C1.18.4) (CL488-65208, Clone: C1.18.4) (blue). Cells were not fixed. Lymphocytes were gated.