

For Research Use Only

CoraLite® Plus 488-conjugated ADK Monoclonal antibody

Catalog Number: CL488-66929



Basic Information

Catalog Number:

CL488-66929

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7466

GenBank Accession Number:

BC003568

GeneID (NCBI):

132

UNIPROT ID:

P55263

Full Name:

adenosine kinase

Calculated MW:

41 kDa

Purification Method:

Protein G purification

CloneNo.:

1G1A5

Recommended Dilutions:

IF-P 1:50-1:500

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF-P

Species Specificity:

Human, Mouse, Rat

Positive Controls:

IF-P : human liver cancer tissue,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human liver cancer tissue using CoraLite® Plus 488 ADK antibody (CL488-66929, Clone: 1G1A5) at dilution of 1:200.