

For Research Use Only

CoraLite® Plus 488-conjugated FTCD Monoclonal antibody



Catalog Number: CL488-66979

Basic Information

Catalog Number: CL488-66979	GenBank Accession Number: BC136395	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 10841	CloneNo.: 1H8A1
Source: Mouse	Full Name: formiminotransferase cyclodeaminase	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG2b	Calculated MW: 541 aa, 59 kDa	Excitation/Emission maxima wavelengths: 488 nm / 515 nm
Immunogen Catalog Number: AG17429	Observed MW: 53-62 kDa	

Applications

Tested Applications:
IF

Species Specificity:
Human, Mouse, Rat, Pig

Positive Controls:

IF : human kidney tissue, human liver cancer tissue

Background Information

Storage


Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human kidney tissue using CoraLite® Plus 488 FTCD antibody (CL488-66979, Clone: 1H8A1) at dilution of 1:200. DAPI (blue).