

For Research Use Only

CoraLite® Plus 488-conjugated PR3/PRTN3 Monoclonal antibody

Catalog Number:CL488-67030



Basic Information

Catalog Number:

CL488-67030

Size:

100ul , Concentration: 1000 ug/ml by

Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG25915

GenBank Accession Number:

BC096183

GeneID (NCBI):

5657

UNIPROT ID:

P24158

Full Name:

proteinase 3

Calculated MW:

256 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

2D1B8

Recommended Dilutions:

IF-P 1:50-1:500

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF-P

Species Specificity:

human

Positive Controls:

IF-P : human tonsillitis tissue,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 PRTN3 antibody (CL488-67030, Clone: 2D1B8) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 PRTN3 antibody (CL488-67030, Clone: 2D1B8) at dilution of 1:200.