

For Research Use Only

CoraLite® Plus 488-conjugated PRTN3 Monoclonal antibody



Catalog Number:CL488-67030

Basic Information

Catalog Number: CL488-67030	GenBank Accession Number: BC096183	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 5657	CloneNo.: 2D1B8
Source: Mouse	Full Name: proteinase 3	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG2b	Calculated MW: 256 aa, 28 kDa	Excitation/Emission maxima wavelengths: 488 nm / 515 nm
Immunogen Catalog Number: AG25915	Observed MW: 28 kDa	

Applications

Tested Applications:
IF
Species Specificity:
Human

Positive Controls:
IF : human tonsillitis tissue,

Background Information

Storage


Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 PRTN3 antibody (CL488-67030, Clone: 2D1B8) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 PRTN3 antibody (CL488-67030, Clone: 2D1B8) at dilution of 1:200.