

For Research Use Only

# CoraLite® Plus 488-conjugated ADRBK1 Monoclonal antibody

Catalog Number:CL488-67067



## Basic Information

Catalog Number:

CL488-67067

Size:

100ul , Concentration: 1000 ug/ml by  
Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG5071

GenBank Accession Number:

BC037963

GeneID (NCBI):

156

UNIPROT ID:

P25098

Full Name:

adrenergic, beta, receptor kinase 1

Calculated MW:

689 aa, 80 kDa

Observed MW:

80 kDa

Purification Method:

Protein A purification

CloneNo.:

2F7A5

Recommended Dilutions:

WB 1:500-1:1000

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

WB

Species Specificity:

Human

Positive Controls:

WB : U-937 cells, THP-1 cells, HL-60 cells

## Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with CL488-67067 (ADRBK1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.