

For Research Use Only

CoraLite® Plus 488-conjugated ADRBK1 Monoclonal antibody



Catalog Number: CL488-67067

Basic Information

Catalog Number: CL488-67067	GenBank Accession Number: BC037963	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 156	CloneNo.: 2F7A5
Source: Mouse	Full Name: adrenergic, beta, receptor kinase 1	Recommended Dilutions: WB 1:500-1:1000
Isotype: IgG2a	Calculated MW: 689 aa, 80 kDa	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG5071	Observed MW: 80 kDa	

Applications

Tested Applications: WB	Positive Controls: WB : U-937 cells, THP-1 cells, HL-60 cells
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with CL488-67067 (ADRB1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.