

For Research Use Only

CoraLite® Plus 488-conjugated GNL3 Monoclonal antibody



Catalog Number: CL488-67169

Basic Information

Catalog Number:

CL488-67169

Size:

100ul

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG7056

GenBank Accession Number:

BC001024

GeneID (NCBI):

26354

UNIPROT ID:

Q9BVP2

Full Name:

guanine nucleotide binding protein-like 3 (nucleolar)

Calculated MW:

62 kDa

Observed MW:

62 kDa

Purification Method:

Protein A purification

CloneNo.:

3B8F10

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data