

For Research Use Only

CoraLite® Plus 488-conjugated AHCYL2 Monoclonal antibody

Catalog Number: CL488-67484



Basic Information

Catalog Number:

CL488-67484

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG19614

GenBank Accession Number:

BC024325

GeneID (NCBI):

23382

UNIPROT ID:

Q96HN2

Full Name:

S-adenosylhomocysteine hydrolase-
like 2

Calculated MW:

611 aa, 67 kDa

Observed MW:

67 and 57 kDa

Purification Method:

Protein A purification

CloneNo.:

1D1F6

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

Human

Positive Controls:

IF/ICC : PC-3 cells,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed PC-3 cells using Coralite® Plus 488 AH CYL2 antibody (CL488-67484, Clone: 1D1F6) at dilution of 1:200. Blue (DAPI).