

For Research Use Only

CoraLite® Plus 488-conjugated AHCYL2 Monoclonal antibody



Catalog Number:CL488-67484

Basic Information

Catalog Number: CL488-67484	GenBank Accession Number: BC024325	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 23382	CloneNo.: 1D1F6
Source: Mouse	Full Name: S-adenosylhomocysteine hydrolase- like 2	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG2a	Calculated MW: 611 aa, 67 kDa	Excitation/Emission maxima wavelengths: 488 nm / 515 nm
Immunogen Catalog Number: AG19614	Observed MW: 67 and 57 kDa	

Applications

Tested Applications: IF	Positive Controls: IF : PC-3 cells,
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed PC-3 cells using CoraLite® Plus 488 AHCYL2 antibody (CL488-67484, Clone: 1D1F6) at dilution of 1:200. Blue (DAPI).