

For Research Use Only

CoraLite® Plus 488-conjugated GSDMB Monoclonal antibody



Catalog Number:CL488-67527

Basic Information

Catalog Number: CL488-67527	GenBank Accession Number: BC025682	Purification Method: Protein G purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 55876	CloneNo.: 2F12C4
Source: Mouse	Full Name: gasdermin B	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG1	Calculated MW: 411 aa, 47 kDa	
Immunogen Catalog Number: AG28332	Observed MW: 45-50 kDa	

Applications

Tested Applications: IF	Positive Controls: IF : HepG2 cells,
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using CoraLite® Plus 488 GSDMB antibody (CL488-67527, Clone: 2F12C4) at dilution of 1:200.