

For Research Use Only

CoraLite® Plus 488-conjugated SNF8 Monoclonal antibody

Catalog Number: CL488-67696



Basic Information

Catalog Number:

CL488-67696

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30384

GenBank Accession Number:

BC008976

GeneID (NCBI):

11267

UNIPROT ID:

Q96H20

Full Name:

SNF8, ESCRT-II complex subunit, homolog (S. cerevisiae)

Calculated MW:

258 aa, 29 kDa

Observed MW:

29 kDa

Purification Method:

Protein G purification

CloneNo.:

2C10F2

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : HeLa cells,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using Coralite® Plus 488 SNF8 antibody (CL488-67696, Clone: 2C10F2) at dilution of 1:200, CL594-Phalloidin (red).