

For Research Use Only

# CoraLite® Plus 488-conjugated RPL6 Monoclonal antibody



Catalog Number:CL488-67729

## Basic Information

<b>Catalog Number:</b> CL488-67729	<b>GenBank Accession Number:</b> BC004138	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 100ul , Concentration: 1000 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 6128	<b>CloneNo.:</b> 3F6F8
<b>Source:</b> Mouse	<b>Full Name:</b> ribosomal protein L6	<b>Recommended Dilutions:</b> IF 1:300-1:1200
<b>Isotype:</b> IgG1	<b>Calculated MW:</b> 33 kDa	<b>Excitation/Emission maxima wavelengths:</b> 488 nm / 515 nm
<b>Immunogen Catalog Number:</b> AG6840	<b>Observed MW:</b> 34 kDa	

## Applications

**Tested Applications:**

IF

**Species Specificity:**

Human, mouse, rat

**Positive Controls:**

IF : HepG2 cells,

## Background Information

Ribosomes are large multiple protein complexes which is composed of more than 80 components and provide a protein-synthesizing machinery for all living cells. Ribosomal protein like 6 (RPL6, also called Taxreb107 for Tax responsive element binding protein 107), one of the ribosomal components, can bind with the Tax responsive element of human T cell leukemia virus (HTLV)-I in addition to be assembled to ribosomes [PMID:12031506]. Overexpressed in gastric cancer could promote cell growth and cell cycle progression at least through upregulating cyclin E expression[PMID:20171175].

## Storage

**Storage:**

Store at -20°C. Avoid exposure to light.

**Storage Buffer:**

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using CoraLite® Plus 488 RPL6 antibody (CL488-67729, Clone: 3F6F8) at dilution of 1:600.