

For Research Use Only

CoraLite® Plus 488-conjugated SCAMP3 Monoclonal antibody

Catalog Number:CL488-67932



Basic Information

Catalog Number: CL488-67932	GenBank Accession Number: BC000161	Purification Method: Protein G purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 10067	CloneNo.: 1C1E6
Source: Mouse	UNIPROT ID: O14828	Recommended Dilutions: IF/ICC 1:50-1:500
Isotype: IgG1	Full Name: secretory carrier membrane protein 3	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG26073	Calculated MW: 38 kDa	
	Observed MW: 33 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : HeLa cells,
Species Specificity: human, mouse, rat	


Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CoraLite® Plus 488 SCAMP3 antibody (CL488-67932, Clone: 1C1E6) at dilution of 1:200.