

For Research Use Only

CoraLite® Plus 488-conjugated SH3BP4 Monoclonal antibody



Catalog Number: **CL488-68045**

Basic Information

Catalog Number: CL488-68045	GenBank Accession Number: BC057396	Purification Method: Protein G purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 23677	CloneNo.: 2G4H3
Source: Mouse	UNIPROT ID: Q9POV3	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Isotype: IgG1	Full Name: SH3-domain binding protein 4	
Immunogen Catalog Number: AG31238	Calculated MW: 963 aa, 108 kDa	
	Observed MW: 110 kDa	

Applications

Tested Applications:
Species Specificity:
Human

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

