

For Research Use Only

# CoraLite® Plus 488-conjugated Acetyl-Histone H3 (Lys23) Recombinant antibody

Catalog Number: **CL488-82838-2**



## Basic Information

Catalog Number:

CL488-82838-2

GenBank Accession Number:

BC066245

Purification Method:

Protein A purification

Size:

100ul , Concentration: 1000 ug/ml by

Nanodrop;

GeneID (NCBI):

8350

CloneNo.:

241137E4

Source:

Rabbit

UNIPROT ID:

P68431

Recommended Dilutions:

IF/ICC 1:50-1:500

Isotype:

IgG

Full Name:

histone cluster 1, H3a

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF/ICC

Positive Controls:

IF/ICC : HeLa cells,

Species Specificity:

human, mouse, rat

## Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 488 Acetyl-Histone H3 (Lys23) antibody (CL488-82838-2, Clone: 241137E4 ) at dilution of 1:200.