

For Research Use Only

CoraLite® Plus 488-conjugated MAZ Recombinant monoclonal antibody

Catalog Number: CL488-82852-3



Basic Information

Catalog Number: CL488-82852-3	GenBank Accession Number: BC041629	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 4150	CloneNo.: 1D5
Source: Rabbit	UNIPROT ID: P56270	Recommended Dilutions: IF/ICC: 1:50-1:500
Isotype: IgG	Full Name: MYC-associated zinc finger protein (purine-binding transcription factor)	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG15200	Calculated MW: 477 aa, 49 kDa	
	Observed MW: 49-51 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : HeLa cells,
Species Specificity: human	

Background Information

Myc-associated zinc-finger protein (MAZ) is a transcription factor with dual roles in transcription initiation and termination. MAZ can activate the expression of a variety of genes such as the c-Myc, insulin, VEGF, and RAS gene family and the serotonin 1A receptor. MAZ plays a critical role in the progression of various cancers including breast, prostate, liver, and so on. Dereglulation of MAZ expression is associated with the progression of pancreatic ductal adenocarcinoma (PDAC)(PMID: 29414775).


Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 488 MAZ antibody (CL488-82852-3, Clone: 1D5) at dilution of 1:200, CL594-phalloidin (red).