

For Research Use Only

CoraLite® Plus 488-conjugated PPFIA3 Recombinant antibody

Catalog Number: CL488-83142-4



Basic Information

Catalog Number:

CL488-83142-4

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1971

GenBank Accession Number:

BC021255

GeneID (NCBI):

8541

UNIPROT ID:

O75145

Full Name:

protein tyrosine phosphatase, receptor type, f polypeptide (PTPRF), interacting protein (liprin), alpha 3

Calculated MW:

1194 aa, 133 kDa

Observed MW:

140 kDa

Purification Method:

Protein A purification

CloneNo.:

230459F6

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human, mouse, rat

Background Information

PPFIA3, also named as KIAA0654, belongs to the liprin family and Liprin-alpha subfamily. PPFIA3 may regulate the disassembly of focal adhesions. It may localize receptor-like tyrosine phosphatases type 2A at specific sites on the plasma membrane, possibly regulating their interaction with the extracellular environment and their association with substrates. PPFIA3 exhibits more than 5 fold accumulations in the monoribosome fraction. PPFIA3 is decreased at the protein level by the miR-99 family but not in the luciferase reporter with the 3'UTR of PPFIA3. (PMID:21212412)

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ HEK-293 cells were intracellularly stained with 0.25 ug CoraLite® Plus 488 Ppfia3 Recombinant Antibody (CL488-83142-4, Clone:230459F6)(red), or 0.25 ug CoraLite® Plus 488 Rabbit IgG Isotype Control Recombinant Antibody (CL488-98136, Clone: 240953C9) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).