

For Research Use Only

CoraLite® Plus 488-conjugated MAPK10 Recombinant antibody

Catalog Number:CL488-84875-5



Basic Information

Catalog Number:

CL488-84875-5

GenBank Accession Number:

BC065516

Purification Method:

Protein A purification

Size:

100ul , Concentration: 1000 ug/ml by

Nanodrop;

GeneID (NCBI):

5602

CloneNo.:

242364A10

Source:

Rabbit

UNIPROT ID:

P53779

Recommended Dilutions:

IF/ICC: 1:50-1:500

Isotype:

IgG

Full Name:

mitogen-activated protein kinase 10

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

Immunogen Catalog Number:

AG11785

Calculated MW:

464 aa, 53 kDa

Applications

Tested Applications:

IF/ICC

Positive Controls:

IF/ICC : HepG2 cells,

Species Specificity:

human

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using CoraLite® Plus 488 MAPK10 antibody (CL488-84875-5, Clone: 242364A10) at dilution of 1:200, CL594-Phalloidin (red).