

For Research Use Only

CoraLite® Plus 488-conjugated LARP2 Recombinant antibody

Catalog Number: CL488-84898-3



Basic Information

Catalog Number:

CL488-84898-3

Size:

100ul, Concentration: 1000 ug/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11043

GenBank Accession Number:

BC030516

GeneID (NCBI):

55132

UNIPROT ID:

Q659C4

Full Name:

La ribonucleoprotein domain family,
member 2

Calculated MW:

201 aa, 24 kDa

Purification Method:

Protein A purification

CloneNo.:

242486E9

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : U2OS cells,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using CoraLite® Plus 488 LARP2 antibody (CL488-84898-3, Clone: 242486E9) at dilution of 1:200, CL594-Phalloidin (red).