

For Research Use Only

CoraLite® Plus 488-conjugated RBPMS Recombinant monoclonal antibody

Catalog Number:CL488-85086-5



Basic Information

Catalog Number:

CL488-85086-5

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7256

GenBank Accession Number:

BC003608

GeneID (NCBI):

11030

UNIPROT ID:

Q93062

Full Name:

RNA binding protein with multiple
splicing

Calculated MW:

22 kDa

Observed MW:

28-30 kDa

Purification Method:

Protein A purification

CloneNo.:

24273087

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse, rat

Positive Controls:

IF/ICC : HepG2 cells,

Background Information

RBPMS is a member of the RNA recognition motif family of RNA-binding proteins. It has a single, putative RNA recognition motif in its N-terminus.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using CoraLite® Plus 488 RBPMS antibody (CL488-85086-5, Clone: 242730B7) at dilution of 1:200, CL594-Phalloidin (red).