

For Research Use Only

CoraLite® Plus 488-conjugated FAM13A Recombinant antibody

Catalog Number:CL488-85092



Basic Information

Catalog Number:

CL488-85092

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014883

GeneID (NCBI):

10144

UNIPROT ID:

O94988

Full Name:

family with sequence similarity 13,
member A

Calculated MW:

117 kDa

Purification Method:

Protein A purification

CloneNo.:

242285A8

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse

Positive Controls:

IF/ICC : A431 cells,

Background Information

FAM13A, also named as FAM13A1 and KIAA0914, belongs to the FAM13 family. It is associated with chronic obstructive pulmonary disease.(PMID: 20173748) This antibody is specific to FAM13A.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using CoraLite® Plus 488 FAM13A antibody (CL488-85092, Clone: 242285A8) at dilution of 1:200.