

For Research Use Only

CoraLite® Plus 488-conjugated THRAP3 Recombinant antibody

Catalog Number: CL488-85093-2



Basic Information

Catalog Number:

CL488-85093-2

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005119

GeneID (NCBI):

9967

UNIPROT ID:

Q9Y2W1

Full Name:

thyroid hormone receptor associated
protein 3

Calculated MW:

109 kDa

Observed MW:

109-130 kDa

Purification Method:

Protein A purification

CloneNo.:

242268D8

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : HepG2 cells,

Background Information

THRAP3, also named as TRAP150, plays a role in transcriptional coactivation.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using CoraLite® Plus 488 THRAP3 antibody (CL488-85093-2, Clone: 242268D8) at dilution of 1:200, CL594-phalloidin (red).