

For Research Use Only

CoraLite® Plus 488-conjugated ETS2 Recombinant antibody

Catalog Number: CL488-85113-4



Basic Information

Catalog Number:

CL488-85113-4

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2929

GenBank Accession Number:

BC017040

GeneID (NCBI):

2114

UNIPROT ID:

P15036

Full Name:

v-ets erythroblastosis virus E26
oncogene homolog 2 (avian)

Calculated MW:

469 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

242817H8

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : MCF-7 cells, HeLa cells

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 488 ETS2 antibody (CL488-85113-4, Clone: 242817H8) at dilution of 1:200, CL594-Phalloidin (red).



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using CoraLite® Plus 488 ETS2 antibody (CL488-85113-4, Clone: 242817H8) at dilution of 1:200, CL594-Phalloidin (red).