

For Research Use Only

CoraLite® Plus 488-conjugated Parvalbumin Recombinant monoclonal antibody

Catalog Number: CL488-85819-4



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
CL488-85819-4	NM_002854	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
100ul , Concentration: 1000 ug/ml by Nanodrop;	5816	250158B6
Source:	UNIPROT ID:	Recommended Dilutions:
Rabbit	P20472	IF-P: 1:50-1:500
Isotype:	Full Name:	Excitation/Emission maxima wavelengths:
IgG	parvalbumin	493 nm / 522 nm
Immunogen Catalog Number:	Calculated MW:	
AG30189	12 kDa	

Applications

Tested Applications:	Positive Controls:
IF-P	IF-P : mouse cerebellum tissue,
Species Specificity:	
human, mouse	

Background Information

PVALB is a high affinity calcium ion-binding protein that is structurally and functionally similar to calmodulin and troponin C. PVALB is expressed in high levels only in fast-contracting muscles and at lower levels in brain and several endocrine tissues. It is thought to be involved in muscle relaxation. Monomeric (12 kDa) and dimeric (24 kDa) forms, as well as trimeric Parvalbumin (38 kDa), have been described in a literature (PMID: 19616851).

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse cerebellum tissue using CoraLite® Plus 488 Parvalbumin antibody (CL488-85819-4, Clone: 250158B6) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse cerebellum tissue using CoraLite® Plus 488 Parvalbumin antibody (CL488-85819-4, Clone: 250158B6) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).